

Achieving Lower Detection Limits Easily with the Agilent Multimode Inlet (MMI)

Application Note

All Industries

Authors

Bill Wilson and Chin-Kai Meng
Agilent Technologies, Inc.
2850 Centerville Road
Wilmington, DE 19808
USA

Abstract

This application note discusses three injection techniques: hot splitless, cold splitless, and solvent vent mode available on the Multimode Inlet. The cold splitless and solvent vent mode injections allow analysts to achieve a lower detection limit by making large volume injections (LVI). A total ion chromatogram overlay of 40-ppb pesticide standards from 2- μ L hot splitless, 10- μ L cold splitless and 25- μ L solvent vent illustrates the improvement in signal-to-noise ratios using LVI.



Agilent Technologies

Introduction

A growing number of analysts are exploring large volume injection (LVI) techniques to improve existing analyses. With traditional liquid injection techniques in capillary gas chromatography, most inlets and columns can only handle 1 – 2 μL at a time. Attempts to increase the injection volume can lead to broadened and distorted analyte peaks, large and long solvent peak tails, and saturated or damaged detectors.

The purpose of increasing the injection volume is normally to improve detection limits in trace analysis. By introducing more of the sample to the system, the mass of analyte reaching the detector will be proportionally increased, resulting in larger peak areas and peak heights. If the baseline noise is constant, larger peak heights mean greater signal to noise ratios and lower system detection limits. An additional benefit of LVI is the ability to reduce the amount of sample originally processed. By injecting 10 – 100 times more volume of processed sample and concentrating it in the inlet, the sample preparation can start with 10 – 100 times smaller sample volume and still achieve the same mass of analyte on column. Another advantage of using LVI (solvent vent) is the decrease in solvent that actually reaches the detector. Usually, only 10 – 30% of the injection solvent actually enters the column and makes it to the detector.

LVI can be applied to injection volumes ranging from a few microliters up to 1 mL or more. In most LVI approaches, the sample solvent is evaporated and removed from the inlet system before the analytes are transferred to the separation column. In this way, LVI is similar to nitrogen evaporation or rotary evaporation of the solvent, with the added benefit of being performed in the GC inlet rather than in a fume hood. Analytes that would be lost during nitrogen evaporation may be retained in the inlet and successfully analyzed via LVI. Furthermore, the LVI process can be automated and is reproducible. As in the other evaporation techniques, the LVI approach is a function of the solvent type, the inlet temperature, the vent flow of evaporation gas, and the analyte boiling point. In addition, the inlet pressure during evaporation and the inlet liner have an impact on the rate of solvent removal and analyte recovery. These parameters will be discussed in this application note.

Experimental

MMI Operational Modes

The Agilent Multimode Inlet (MMI) uses the same liners and consumables as a standard split/splitless inlet, making it compatible with existing hot split and splitless methods. Its operational modes include: Hot Split/Splitless (also in pulsed

mode), Cold Split/Splitless (also in pulsed mode), Solvent Vent and Direct mode.

Hot Splitless (for 1 – 3 μL injections)

For most analysts considering LVI, their current methods are using hot splitless injection. This proven and reliable sample introduction technique has worked well for almost 40 years; however, it does present some challenges to the sample integrity and to the method developer. First, the inlet must be hot enough to flash vaporize the solvent and analytes so that the resulting vapor cloud can be transferred to the column. The inlet liner volume must be sufficiently large to contain this vapor cloud. If the liner volume is too small, the vaporized sample can overflow the liner and reach reactive surfaces, leading to analyte loss. In addition, the pressure wave generated by the vaporized sample can push back against the incoming carrier gas and enter sensitive pressure and flow control systems. Using the Agilent pressure/flow calculator [1], a 1- μL injection of acetone into an inlet at 240 °C and 14.5 psig expands to 288 μL of gas. Most inlet liners for standard split/splitless inlets have a nominal volume of 1 mL. An increase of injection volume to only 3.5 μL under these conditions creates a vapor cloud of 1 mL which could easily overflow the inlet liner.

Hot splitless injection also creates a challenging environment for thermally unstable or labile analytes. Compounds such as the organochlorine pesticides DDT and endrin can rearrange to form breakdown compounds. This process is accelerated with the inlet temperatures normally used to analyze them. Effective chemical deactivation of the liner can minimize analyte breakdown. However, high inlet temperatures can decrease the lifetime of deactivated liners.

Another challenge created by hot splitless injection is the opportunity for needle fractionation or analyte discrimination. The needle temperature increases as the sample is being transferred from the syringe to the inlet because the needle is in contact with the septum. The rise in needle temperature can cause the solvent to "boil" away and deposit high boiling analytes inside the needle. To avoid this fractionation problem, some analysts load a solvent plug into the syringe first and then draw up the desired sample volume (available in 7693A Automatic Liquid Sampler). The thought is that the solvent plug will wash any deposits into the inlet. An effective way to address this problem is to make a high speed injection. This minimizes the time the needle is in contact with the septum and the time the sample touches the needle. Even with these issues, hot splitless injection is a well-accepted technique. An alternative technique, such as cold splitless can address these concerns and improve the analysis results.

Cold Splitless (for 1 – 10 µL injections)

MMI's versatile temperature programmability allows it to perform cold split and splitless analyses. In cold splitless mode, the MMI is cooled to a temperature below the normal boiling point of the sample solvent so that when the sample is injected, no vaporization takes place. The injection is simply a liquid transfer from the syringe to the inlet. Once the syringe is removed from the inlet, the inlet is heated to vaporize the sample and transfer it to the column. The solvent vaporizes first and moves to column, allowing analyte focusing to take place as in normal hot splitless injections. The analytes subsequently vaporize and move to the column. The main advantage is that the analytes vaporize at the lowest possible inlet temperature, rather than at a constant high temperature. This minimizes thermal degradation while still allowing a wide range of analytes to vaporize. Cold splitless operations also do not thermally stress the liner as harshly as hot splitless does, prolonging its usable life. Cold splitless can also extend the amount of sample that can be injected in some cases. If a slow inlet temperature program is used, the solvent can be vaporized slowly and will not overflow the liner volume. As long as the analytes can be refocused on the column, slow inlet temperature programs cause no detrimental effects to the chromatography.

Solvent Vent (for 5 – 1000 µL injections)

The solvent vent mode is the method which enables MMI to do LVI of more than 5 µL. In solvent vent mode, the inlet is kept at a low initial temperature during sample injection. Pneumatically, the inlet is in split mode with a low inlet pressure. The flow of gas through the inlet liner and out to vent removes the evaporating solvent. The sample is injected slowly so that the incoming liquid is deposited on the liner wall and the solvent evaporates at a similar rate. Once the entire sample has been injected, the inlet switches to a splitless mode for analyte transfer. The inlet is then heated to vaporize the concentrated sample and any remaining solvent and the vapor is transferred to the column. After a sufficient period to ensure the sample transfer, the inlet is then switched to a purge mode to allow any remaining material in the inlet liner to be vented. During the sample injection and solvent venting period, the GC oven has been held at an appropriate temperature to allow the solvent to refocus the analytes on the column. When this refocusing is complete, the oven is then programmed to perform the separation.

LVI Method Development

An effective procedure for developing an LVI method on a MMI is to run the existing method first to determine peak areas for a small volume injection. Such results serve as a baseline for evaluating the LVI method performance. The next step is to switch to the solvent vent mode with a slightly larger injection volume (for example, 2 to 5 times larger). By comparing the resulting peak areas and accounting for the increased injection volume, the analyte recovery can be calculated and conditions can be further optimized.

Backflush

A traditional bakeout step for removing late eluters can be very time consuming for samples with complicated matrices, even as long as the analysis time. Capillary flow devices (in this case, a purged ultimate union) provide backflush [2, 3] capability. "Backflush" is a term used for the reversal of flow through a column such that sample components in the column are forced back out the inlet end of the column. By reversing column flow immediately after the last compound of interest has eluted, the long bake-out time for highly retained components can be eliminated. Therefore, the column bleed and ghost peaks are minimized, the column will last longer, and the MS ion source will require less frequent cleaning. The split vent trap may require replacement more frequently than usual.

Instrument Parameters

GC	Agilent 7890A
MS	Agilent 5975C MSD
Column	HP-5MS UI, 15 m × 0.25 mm × 0.25 µm (19091S-431UI), from inlet to purged union
MMI	Constant pressure (~18 psi), chlorpyrifos-methyl RT locked to 8.297 min, 2 psi at post run for backflush
MMI liner	Double taper deactivated, Helix (5188-5398)
Septum purge	3 mL/min
Purged Union	4 psi; 70 psi at post run for backflush
Restrictor	0.7 m × 0.15 mm deactivated fused silica tubing (from purged union to MSD)
Syringes	10 µL, for splitless injections (5181-3354) 50 µL, for solvent vent mode (5183-0318)
ALS	Agilent 7693A
MS parameters	
Solvent delay	2.5 min
Gain factor	1
Mass range	44–550
Threshold	0
Samples	2
Tune file	atune.u

Oven

Initial temperature	70 °C
Initial hold time	1 min
Rate 1	50 °C/min
Temperature 1	150 °C
Hold time	0 min
Rate 2	6 °C/min
Temperature 2	200 °C
Hold time	0 min
Rate 3	16 °C/min
Temperature 3	280 °C
Hold time	5 min
Total runtime	20.933 min
Post run	5 min (for backflush)
Oven post run temp	280 °C

Sample: 40-ppb pesticide standards in acetone (for a list of compounds, see Figure 5).

Multimode Inlet (MMI)

Parameter	Hot Splitless	Cold Splitless	Solvent Vent
Initial temperature	280 °C	30 °C	35 °C
Initial time	–	0.01 min	0.35 min
Rate 1	–	700 C/min	700 °C/min
Final temperature	–	320 °C	320 °C
Vent flow	–	–	150 mL/min
Vent pressure	–	–	5 psig
Vent time	–	–	0.33 min (from calculator, Figure 3)
Purge time	0.75 min	1.25 min	1.5 min
Purge flow	50 mL/min	50 mL/min	50 mL/min
Injection volume	2 µL	10 µL	25 µL
Injection speed	Fast	Fast	75 µL/min (from calculator, Figure 3)
Cryo	–	On (liquid CO ₂)	On (liquid CO ₂)
Cryo fault detection	–	On	On
Cryo use temperature	–	125 °C	125 °C
Time out detection	–	On (15 min)	On (15 min)

The parameters for the 25-µL Solvent Vent injection were determined with the Solvent Elimination Calculator integrated in the ChemStation. This calculator was designed to help determine reasonable starting conditions for LVI methods. When the MMI is put into the PTV Solvent Vent mode, an additional button appears in the inlet screen, shown in Figure 1.

In the first screen of the Solvent Elimination Calculator (Figure 2), the sample solvent and desired injection volume are selected and entered. The calculator "knows" the syringe currently installed and will only allow 50% of that volume to be injected at once. Larger injection volumes can be entered into the calculator but the injection volume will not be downloadable. The calculator also requests the boiling point of the earliest eluting analyte, as this allows the initial inlet temperature to be selected. If the boiling point is unknown, the temperature should be left at 150 °C as this will work for a wide range of analytes.




Figure 1. Multimode Inlet "Solvent Elimination Calculator" imbedded in ChemStation for easy method development.




Figure 2. Select solvent of choice and enter the injection volume to start the calculation.

Figure 3 shows the calculation screen. The calculator uses an initial set of inlet conditions to determine the solvent elimination rate according to fundamental theory [4]. This "Elimination Rate" does not account for other factors (for example, local cooling due to solvent evaporation) specific to LVI and is normally faster than that determined from practical experience. The "Suggested Injection Rate" does consider these factors and is designed to leave a small amount of solvent in the liner at the end of the venting period. This solvent serves as a liquid "trap" for the more volatile analytes and promotes their recovery. The "Suggested Vent Time" is determined by dividing the injection volume by the "Suggested Injection Rate."

Several variables for determining elimination rate can be set by the user in the lower portion of the window. A small change in inlet temperature has a significant impact on elimination rate. Vent flow has a linear effect such that a decrease by a factor of two in vent flow gives an equal decrease in elimination rate. As the vent pressure decreases, the elimination rate increases. Bear in mind that the vent pressure also impacts the amount of solvent that reaches the column during venting. As the vent pressure is increased, more solvent is loaded onto the column before the analytes are transferred. Finally, the type of solvent, specifically its normal boiling point, has a substantial impact on the elimination rate.

Figure 3. The calculator calculates the injection rate and vent time according to the selected inlet temperature and vent flow.

The download screen in Figure 4 shows all of the method changes that are downloaded to the edit parameters screen. The check boxes allow the user to accept (by checking) or reject any of these parameters. The oven initial temperature and hold times are not automatically checked in case the current method requires these values to be unchanged (for example, a Retention Time Locked method).




Figure 4. Confirm values suggested by the Calculator and download to ChemStation.

Results and Discussion

Figure 5 compares the responses of a 40-ppb standard solution from three injection modes.

The bottom total ion chromatogram (TIC) is a typical 2- μ L hot splitless injection. Some of the 40-ppb pesticides are barely visible (80 pg each on column). The middle TIC is from a 10- μ L cold splitless injection. The MMI starting temperature was

30 °C. In this TIC, the on column amount for each analyte is 400 pg. Lastly, the top TIC is from a 25- μ L solvent vent injection with MMI starting temperature at 35 °C. In this TIC, the signal-to-noise ratio is significantly better than the TIC from hot splitless injection (bottom TIC), as noted in the Introduction section. The peak shape and resolution are maintained, even with the 25- μ L injection volume. This implies that the solvent was mostly eliminated during the injection.




Figure 5. Overlay of total ion chromatograms (TICs) from three injection modes, plotted on the same scale.

Conclusion

The new Agilent Multimode Inlet (MMI) has the same form factor and uses the same consumables (for example, liners, o-rings and septa) as the existing split/splitless inlet, allowing existing hot splitless methods to be replicated. In addition, the temperature programmability permits both cold splitless and large volume injection (LVI) methods for improved detection limits. An integrated Solvent Elimination Calculator provides a complete set of initial conditions for easy LVI method development. The application results show a significant signal-to-noise improvement (lower detection limits) comparing the 25- μ L solvent vent injection to the 2- μ L hot splitless injection.

References

1. Agilent Pressure/Flow Calculator Included in the Instrument Utility DVD, available with each gas chromatograph and MMI accessory kit.
2. Chin-Kai Meng, "Improving Productivity and Extending Column Life with Backflush," Agilent Technologies publication, 5989-6018EN, December 2006.
3. Matthew Klee, "Simplified Backflush Using Agilent 6890 GC Post Run Command," Agilent Technologies publication, 5989-5111EN, June 2006.
4. J. Stanieski and J. Rijks, *Journal of Chromatography* 623 (1992) 105-113.

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2009
Printed in the USA
June 18, 2009
5990-4169EN



Agilent Technologies